



Thermo Fisher launches collection kit to streamline saliva testing

30 June 2021 | News

The individually packaged and barcoded kits require just 1 mL of saliva

To support large-scale research and surveillance programs, Thermo Fisher Scientific has introduced the Thermo Scientific SpecIMAX Saliva Collection Kit.

While most approved methods for SARS-CoV-2 saliva collection use specialized tubes and buffers for RNA stabilization and virus inactivation, researchers have found that this expensive approach may not always be required.

The SpecIMAX Saliva Collection Kit is designed to provide a widely available, cost-effective solution for raw saliva collection to simplify saliva-based testing, which is necessary for SARS-CoV-2 mapping and epidemiological studies on a global scale.

The easy-to-use self-collection kits facilitate clean saliva transfers, reduce risk of cross-contamination when paired with liquid handlers, and fit seamlessly into viral RNA extraction and direct-to-PCR downstream automation workflows for high throughput surveillance testing.

The individually packaged and barcoded kits require just 1 mL of saliva and include a convenient funnel and cap to secure the sample for processing. The standardized collection tube size can hold 6 mL sample volume yet requires less refrigeration, incubation, and storage space than other collection kits available on the market.